

# LG2000 Anorectal Treating Device

LG2000 Anorectal Treating Device Series belong to a new generation's therapeutic device for proctologic diseases. Its advantages are: with multiple functions; of high technology; Applied in a wide range; of reliable security with little suffering and having fine treatment effect. The particular patent technology of biological electrical impedance auto-measurement (BEIM) successfully settles medical difficulty of haemorrhage after operation by other treatment methods.

## User-friendly touch screen interface



## Scope of application

### Department of coloproctology:

Internal, external and mixed hemorrhoids, anal fissure, anal fistula, proctopolypus, anal abscess, hypertrophy of anal papilla, condyloma acuminatum, etc.

## Main clinical characteristics

- Rapid curative treatment, each hemorrhoid only needs 3-5 seconds;
- No haemorrhage after surgery with little suffering and fine curative effect;
- Make tissues dry and hard without carbonization. Bound of dry and hard part is clear with high security;
- No need to change dressing, no hospitalization, no infection, no sequela and complicating disease;

## Treating System Accessories



### BEIM Pencil

Used for excising the uncalled-for hemorrhoid tissues, assisting the electronic pincers to finish the hemorrhoid treatment. It has 3 different tips (blade, needle, ball) and 9 adjustable positions for output power to meet different demands. Display pencil working power in 0, 5, 15, 20, 25, 30, 35, 40W, total 9 grades.



### BEIM Pincers

BEIM Pincers intelligently measures hemorrhoids, and then make the base of the hemorrhoid dry and hard without carbonization. It acts as sealing machines to seal off the base, no blood can flow out.



### BEIM Forceps

Arresting bleeding instantly for larger blood vessels, making patients timely free from trouble after treatment.

## Configuration

- Main unit (1 set)
- Treating accessories
  1. High-frequency BEIM pincer (2 pcs)
  2. Pincer cable (1 pc)
  3. High-frequency Foot-control ESU Pencil (1 pc)
  4. HF Pencil electrodes (2 pcs)
  5. High-frequency Bipolar Forceps (1 pc)
  6. Forceps cable (1 pc)
  7. Foot switch (1 pc)

## Disposable use Bipolar pincer/forceps

- EO sterilization ,clean and sanitary
- Disposable use,no cross-infection risk
- Increase the profit of hospital or clinic

